

Tool Mate® Fluids

- Fluids & Equipment -

TECHNICAL DATA SHEET

MODEL MP-2

METALLOGRAPHIC SPECIMEN GRINDING & POLISHING MACHINE



Model MP-2 Metallographic specimen grinding & polishing machine is equipped with double discs, and it can finish grinding or polishing with two different speeds. The rotating speed of left disc (pre-grinding disc) is 450 rpm and right disc (polishing disc) is 600 rpm. With double disc and different grinding and polishing materials, it can produce rough grinding, fine grinding, rough polishing and finishing polishing. This machine features easy operation and super value. It is the ideal specimen preparing instrument for use in factories, scientific research institutes and laboratories.

Technical specifications

Diameter of grinding disc: Ø230mm, rotating speed 450 rpm (or customized)

Diameter of polishing disc: Ø200mm, rotating speed 600 rpm (or customized)

Motor power: 370W

Power supply: Single phase, 220V, 50Hz

Dimensions: 70 × 60 × 28 cm Net weight: 40 Kg

For technical inquiries, please email info@toolmatefluids.com

TOOL MATE®
Fluids

1A, Level 2, 802 Pacific Highway
PO Box Locked Bag 1009
Gordon, NSW Australia 2072
PH: +61 2 9844 5457 Fax: +61 2 9844 5445
WEB: www.metalworkingfluids.com.au

